



Energy flow

www.lenyolife.com



HOTEL SOLUTIONS

Attuned to the guests





“I travel a great deal both as a company owner and instructor. Although I am invited to the best hotels, I find it difficult to acclimatize to strange surroundings, which includes adjusting my nervous and digestive system after a long journey. I can also tell by the look of my colleagues and other invited participants just how difficult it is to stay fresh in body and soul, and how hard it is to get some refreshing sleep in a new environment.

This inspired me to develop program packages that can also be applied in similar situations.

It is a great challenge to create a homelike feel in a strange place.”


Gábor Lednyiczky

owner of Hippocampus-Brt Kft, producer of Lenyo devices,
Developer of Bioregulation Therapy (BRT)

The relationship of homeodynamics and BRT

Forces of nature:

- electro-magnetism
- strong interaction
- weak interaction
- gravity



What is homeodynamics? It is the adaptational ability of living organisms to changing external and internal conditions, which ensures their own relative biological stability.

The body exchanges matter and energy with its environment. These “ingredients” are provided for the living organism by self-supporting functions, which are collectively called **self-regulating processes**.

The system of homeodynamics is controlled by electro-magnetic waves

The scientists featured in the scientific history section of our website have proven that the system of homeodynamics is controlled by electro-magnetic waves. The external and internal electro-magnetic connections of the body are collectively called bio-communication.

Bio-regulation (BRT) Technology

In the age of quantum physics, it has been proven by scientific research that the functioning and co-ordination of the cells, tissues, and organs that make up the human body are based on electro-magnetic communication. On average, 7000 chemical reactions take place in our cells in every second. The self-regulation of these processes controlled by electro-magnetic principles is called bio-regulation. BRT is a method that uses pulsed electro-magnetic signals to strengthen the control signals of homeodynamics.

REFERENCES



In the Cancun hotel zone, the B-Well clinic was the first to implement the LENYO Sanocenter solution

Why choose LENYO?

It is a new wellness service which provides an advantage on the market.



Several satisfied, returning guests.



**Easy to use,
economic to operate.
Quick return on investment.**





Extra source of income (dependent on the solution, either included in the room price or offered as a wellness service).



Can be complimentary with any other treatment.



Ongoing technical support. Protocols provided online, tailored to the needs of the hotel.



Flexibility to satisfy the needs of guests for all ages.

“I have never slept so well anywhere else”

According to studies made across several countries, almost every other person struggles with some degree of sleep problem.



“Jet lag is just a memory for me.

I work as a pilot. Even after a long and tiring flight, I adjust quickly to local conditions. I have bought a LENYO device for myself. It would be great if this service was also available at hotels.”

Tibor Vigh, pilot





“I feel fantastic!

Previously, during vacation, I always felt unwell. I’ve struggled with symptoms such as acute exhaustion and irritability... I couldn’t even enjoy the excellent services of the highest quality hotels. Thanks to the LENYO solution, now I enjoy every minute of my vacations.”

“My digestion works well even in a strange environment.

I am an athlete participating in world class sports events. It is important for me to perform well anywhere in the world. Sometimes digestion problems held me back.

That is a thing of the past now.

I have chosen LENYO.”

Erika Csomor, Hungarian long distance runner, duathlon, triathlon and multiple times Ironman



“Active wellness? I enjoy every moment of it.

During vacations, I am often tired, my muscles are sore, and my back and knees ache from the unusual amount of activity (swimming, sightseeing). In this hotel I feel like a new person. I wake up easily and feeling refreshed. I use the LENYO service regularly.”



“I find it easy to focus during long meetings and presentations.

The secret to a successful contract is to be fresh and well rested at meetings. I owe a lot to the hotel. I find it easy to focus at conferences; I hold my presentations feeling fresh and confident. The LENYO service proved to be a good investment.”

LENYO's room devices offer ideal solutions to accommodation establishments of any size, from small boarding houses to the largest hotels. The service may be installed in rooms, included in the room prices, or offered as a separate service (rental). The bed-sized and approximately 5 cm thick signal transmitting mattress can be placed on top of the standard mattress, and easily secured using a mattress protector.

LENYO Sleep – Because comfort during sleep is important.

Easy to use

Guests can easily start the device by pressing a button. It contains two programs: one short (1 hour 6 minutes and 17 seconds) and one long (8 hours 23 minutes and 15 seconds). Guests can use the short program for siesta during the day and the long program during the night. The device may be used with two antennas, for the same program, simultaneously.




LENYO Fractal – Simple yet excellent

Simple

The program packages and signal intensity are easy to select using the buttons on the front panel of the device. It may be used with two easily connectible signal transmitting mattresses simultaneously. This way, a more intense effect is attainable for a single person; alternatively, two people may use the device at the same time, with the same programme package, making it cost efficient.

Excellent

LENYO Fractal's 32 (10-20 minutes) program packages cover the most common and frequent problems that cause unwellness.



Lenyo SanoCenter (LSC) – unique on a global scale

Centralized controls

It is a unique system, controllable from a central computer, for wellness and spa hotels, where it is important to serve the maximum number of guests simultaneously, with a minimum number of staff. This customized service may be distributed from the front desk to a maximum of 64 rooms, using a single controller. It stores previously used protocols as needed, and program combinations may be set as much as 2 weeks earlier, making them easy to offer in packages. The system is flexible to suit the hotel owner's needs.

Many ways of selling

It may be included in the room prices or sold as a service. Guests can select the best program from a list according to their needs (e.g.: sleep, detoxification, digestion, joints, refreshing, concentration, relaxation, weather sensitivity...) in their own rooms, then contact the front desk.

The solution may also be applied in the relaxation and conference facilities of the hotel. The signal transmitters can be built into the relaxation chairs.

A wide range of programs

160 normal and 40 short program packages are included. Their duration varies between 10 and 72 minutes. From these, protocols (combinations of program packages) can be created and stored, as required, for improving the guests' well-being.



The wellness devices recommended by LENYO for wellness centres can be controlled by computer. These are ideal for use in wellness and spa hotels, as a special service of the wellness department, especially where guests are taken care of individually.

LENYO Lux Plus – multi-faceted and programmable

Programmable

Lenyo Lux Plus (LLP) is a wellness device that can be controlled by computer. It stores previously used protocols for an unlimited number of users, which both eases tracking and simplifies use. Preset protocols can be controlled with a single button.

Multi-faceted

96 normal and 40 short program packages are included, their duration varying between 10 and 72 minutes. From these, protocols (combinations of program packages) can be created and stored for potential problems. Two users may run the same program package simultaneously; or, for achieving a more intense effect, a single user may connect 3 antennas at the same time.



LENYO CellCom – Attuned to the individual!

High-tech


The most advanced of the LENYO devices is LENYO CellCom (LCC), which is an endogenous device that can be controlled by computer. Taking into account the changes of the individual's body, it selectively controls the electro-magnetic field in the 1Hz – 1 MHz range. It has more than 500 automatic program packages.

Bio-feedback

Using LCC, a frequency-sensitivity map can be generated for the individual, in the 1 Hz – 1MHz frequency range. The user can observe sensitivity when using special frequency bands through changes in the pulse.

Individual application

Dependent on purpose, protocols (combinations of program packages) suited to the individual can be created and stored.



The devices are CE-marked wellness devices. Their beneficial health effects are not proven.



Energy flow

Hippocampus-Brt Kft.

R&D and Sales Office:
1082 Budapest, Hajós u. 32, Hungary
E-mail: lenyo@lenyolife.com
Tel: + 36 70 3827291 or + 36 70 335 0586
www.lenyolife.com

Distributing partner:

